

TCEB2 Rabbit mAb

Catalog No: #60257

Package Size: #60257-1 50ul #60257-2 100ul

Description

Product Name	TCEB2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	TCEB2 Antibody detects endogenous levels of total TCEB2
Immunogen Description	A synthesized peptide derived from human TCEB2
Other Names	EloB; Elongin B; SIII; SIII p18; Tceb2;
Accession No.	Uniprot:Q15370
Calculated MW	16kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Background

SIII, also known as elongin, is a general transcription elongation factor that increases the RNA polymerase II transcription elongation past template-encoded arresting sites.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.