

RNF7 Rabbit mAb

Catalog No: #60261

Package Size: #60261-1 50ul #60261-2 100ul


Description

Product Name	RNF7 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	RNF7 Antibody detects endogenous levels of total RNF7
Immunogen Description	A synthesized peptide derived from human RNF7
Other Names	CKBBP1; Rbx2; RNF7; ROC2; SAG;
Accession No.	Uniprot:Q9UBF6
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

Images



Western blot analysis of RNF7 expression in HepG2 cell lysate.

Background

Probable component of the SCF (SKP1-CUL1-F-box protein) E3 ubiquitin ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins involved in cell cycle progression, signal transduction and transcription.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.