

BIN1 Rabbit mAb

Catalog No: #60262

Package Size: #60262-1 50ul #60262-2 100ul

Description

Product Name	BIN1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse
Specificity	BIN1 Antibody detects endogenous levels of total BIN1
Immunogen Description	A synthesized peptide derived from human BIN1
Other Names	AMPH2; Amphiphysin 2; Amphiphysin II; Amphiphysin like protein; AMPHL; Bin1; Box Dependand MYC Interacting Protein 1; Bridging integrator 1; Myc box dependent interacting protein 1; SH3P9;
Accession No.	Uniprot:O00499
Calculated MW	50,60,65kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

Background

May be involved in regulation of synaptic vesicle endocytosis. May act as a tumor suppressor and inhibits malignant cell transformation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.