

NDUFAB1 Rabbit mAb

Catalog No: #60264

Package Size: #60264-1 50ul #60264-2 100ul


Description

Product Name	NDUFAB1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP FC
Species Reactivity	Human
Specificity	NDUFAB1 Antibody detects endogenous levels of total NDUFAB1
Immunogen Description	A synthesized peptide derived from human NDUFAB1
Other Names	ACP; FASN2A; ndufab1; SDAP;
Accession No.	Uniprot:O14561
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 IP 1:50 FC 1:50

Images



Western blot analysis of NDUFAB1 expression in A431 cell lysate.

Background

Carrier of the growing fatty acid chain in fatty acid biosynthesis in mitochondria. Accessory and non-catalytic subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I), which functions in the transfer of electrons from NADH to the respiratory chain.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.