

CD275 Rabbit mAb

Catalog No: #60271

Package Size: #60271-1 50ul #60271-2 100ul


Description

Product Name	CD275 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	CD275 Antibody detects endogenous levels of total CD275
Immunogen Description	A synthesized peptide derived from human CD275
Other Names	B7H2; B7RP1; CD275; GL50; ICOSL; Icoslg; LICOS;
Accession No.	Uniprot:O75144
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of CD275 in HUVEC cell lysate.

Background

Ligand for the T-cell-specific cell surface receptor ICOS. Acts as a costimulatory signal for T-cell proliferation and cytokine secretion; induces also B-cell proliferation and differentiation into plasma cells.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.