

RalA Rabbit mAb

Catalog No: #60272

Package Size: #60272-1 50ul #60272-2 100ul


Description

Product Name	RalA Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	RalA Antibody detects endogenous levels of total RalA
Immunogen Description	A synthesized peptide derived from human RalA
Other Names	Ral a; RAL; RALA;
Accession No.	Uniprot:P11233
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Images



Western blot analysis of RalA in MCF7 cell lysate.

Background

Multifunctional GTPase involved in a variety of cellular processes including gene expression, cell migration, cell proliferation, oncogenic transformation and membrane trafficking. Accomplishes its multiple functions by interacting with distinct downstream effectors.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.