

Cathepsin H Rabbit mAb

Catalog No: #60273

Package Size: #60273-1 50ul #60273-2 100ul


Description

Product Name	Cathepsin H Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human
Specificity	Cathepsin H Antibody detects endogenous levels of total Cathepsin H
Immunogen Description	A synthesized peptide derived from human Cathepsin H
Other Names	ACC4; ACC5; cathepsin B3; cathepsin BA; CPSB; CTSH;
Accession No.	Uniprot:P09668
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

Images



Western blot analysis of Cathepsin H in HepG2 cell lysate.

Background

Important for the overall degradation of proteins in lysosomes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.