

Ube1L / UBA7 Rabbit mAb

Catalog No: #60275

Package Size: #60275-1 50ul #60275-2 100ul


Description

Product Name	Ube1L / UBA7 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC FC
Species Reactivity	Human Mouse Rat
Specificity	Ube1L / UBA7 Antibody detects endogenous levels of total Ube1L / UBA7
Immunogen Description	A synthesized peptide derived from human Ube1L / UBA7
Other Names	UBA1B; UBA7; UBE1L; UBE2; Ubiquitin activating enzyme 2; Ubiquitin-activating enzyme 7; Ubiquitin-activating enzyme E1 homolog;
Accession No.	Uniprot:P41226
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 FC 1:50

Images



Western blot analysis of Ube1L / UBA7 in HepG2 cell lysate.

Background

Activates ubiquitin by first adenylating with ATP its C-terminal glycine residue and thereafter linking this residue to the side chain of a cysteine residue in E1, yielding an ubiquitin-E1 thioester and free AMP.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.