

SDHC Antibody

Catalog No: #60284

Package Size: #60284-1 50ul #60284-2 100ul


Description

Product Name	SDHC Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthesized peptide derived from human SDHC
Conjugates	Unconjugated
Target Name	SDHC
Other Names	CYBL; PGL3; QPs1; SDH3; sdhC;
Calculated MW	17kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Images



All lanes use the Antibody at 1:6K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:6K dilution for 1 hour at room temperature.

Background

Membrane-anchoring subunit of succinate dehydrogenase (SDH) that is involved in complex II of the mitochondrial electron transport chain and is responsible for transferring electrons from succinate to ubiquinone (coenzyme Q).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.