

ACADVL/VLCAD Rabbit mAb

Catalog No: #60287

Package Size: #60287-1 50ul #60287-2 100ul

Description

Product Name	ACADVL/VLCAD Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB,IHC,IP
Species Reactivity	Human
Immunogen Description	A synthesized peptide derived from human ACADVL/VLCAD
Conjugates	Unconjugated
Calculated MW	65kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 IP 1:50

Background

Active toward esters of long-chain and very long chain fatty acids such as palmitoyl-CoA, myristoyl-CoA and stearoyl-CoA. Can accommodate substrate acyl chain lengths as long as 24 carbons, but shows little activity for substrates of less than 12 carbons.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.