

Calcitonin mouse mAb

Catalog No: #61044

Package Size: #61044-1 50ul #61044-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Calcitonin mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Clone No.	MR4541
Isotype	IgG
Applications	IHC,IF,FC
Species Reactivity	Human
Conjugates	Unconjugated
Storage	Store at -20 Λ C

Application Details

IHC 1:100-1:200 IF 1:100-1:200

Images




Immunocytochemistry/ Immunofluorescence Calcitonin antibody (61044)

ICC/IF staining of Calcitonin antibody in TT cells. Cells were fixed with 4% Paraformaldehyde permeabilized with 0.1% Triton X-100.

Samples were incubated with 61044 at a working dilution of 1/100. The secondary antibody was Alexa FluorB 647 goat anti mouse, used at a dilution of 1/500.

Nuclei were counterstained with DAPI.



Formalin-fixed, paraffin-embedded human thyroid medullary carcinoma tissue stained for Calcitonin using 61044 at 1/100 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.