

Anti-VKORC1 Antibody

Catalog No: #62003

Package Size: #62003-1 50ul #62003-2 100ul

Description

Product Name	Anti-VKORC1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Human
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the N-terminal region of human VKORC1. The exact sequence is proprietary.
Conjugates	Unconjugated
Calculated MW	Predicted band size: 18 kD; Observed band size: 15 kD
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide

Application Details

WB 1 μ O 500 - 1 μ O 1000

Images



Western blot analysis of VKORC1 expression in human lung
(A)
whole cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.