

B-Myb Polyclonal Antibody

Catalog No: #62006

Package Size: #62006-1 50ul #62006-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	B-Myb Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IP;IF;ELISA
Species Reactivity	Human;Mouse;Rat;Rice
Specificity	B-Myb Polyclonal Antibody detects endogenous levels of B-Myb protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human B-Myb.
Conjugates	Unconjugated
Target Name	MYBL2
Other Names	MYBL2;BMYB;Myb-related protein B;B-Myb;Myb-like protein 2
Calculated MW	80kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500 - 1:2000. IF 1:200 - 1:1000. ELISA: 1:20000

Images



Western blot analysis of lysates from A549 cells, using B-Myb Antibody. The lane on the right is blocked with the synthesized peptide.



Immunofluorescence analysis of HeLa cells, using B-Myb Antibody. The picture on the right is blocked with the synthesized peptide.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.