

IRG1 Rabbit Monoclonal Antibody

Catalog No: #62023

Package Size: #62023-1 50ul #62023-2 100ul

Description

Product Name	IRG1 Rabbit Monoclonal Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purification	Affinity Purification
Applications	WB,IP
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant protein of human IRG1
Conjugates	Unconjugated
Other Names	CAD; ACOD1
Calculated MW	53kDa
Formulation	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Storage	Store at -20°C/1 year

Application Details

WB 1:1000, IP 1:20

Note: This product is for in vitro research use only and is not intended for use in humans or animals.