

EOMES Polyclonal Antibody

Catalog No: #62031

Package Size: #62031-1 50ul #62031-2 100ul



Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	EOMES Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;ELISA
Species Reactivity	Human;Mouse
Specificity	EOMES Polyclonal Antibody detects endogenous levels of EOMES protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human EOMES.
Conjugates	Unconjugated
Target Name	EOMES
Other Names	EOMES;TBR2;Eomesodermin homolog;T-box brain protein 2;T-brain-2;TBR-2
Calculated MW	73kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500 - 1:2000. ELISA: 1:10000.

Images



Western blot analysis of lysate from HT29 cells, using EOMES antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.