

SEMA3A Polyclonal Antibody

Catalog No: #62038

Package Size: #62038-1 50ul #62038-2 100ul


Description

Product Name	SEMA3A Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	SEMA3A Polyclonal Antibody detects endogenous levels of SEMA3A protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human SEMA3A.
Conjugates	Unconjugated
Target Name	SEMA3A
Other Names	SEMA3A;SEMAD;Semaphorin-3A;Semaphorin III;Sema III
Calculated MW	89kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

IHC 1:100 - 1:300. ELISA: 1:40000.. IF 1:50-200

Images



Immunohistochemistry analysis of SEMA3A antibody in paraffin-embedded human tonsil tissue.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.