

ITGB5 Mouse Monoclonal Antibody

Catalog No: #62041

Package Size: #62041-1 50ul #62041-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	ITGB5 Mouse Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	IF,FCM
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human ITGB5 (AA: 50-319) expressed in E. Coli.
Conjugates	Unconjugated
Calculated MW	88kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details


IF 1/50 - 1/200

FCM 1/200 - 1/400

Images



Immunofluorescence analysis of Hela cells using ITGB5 mouse mAb (green).



Flow cytometric analysis of Jurkat cells using ITGB5 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.