

TNNT1 rabbit polyclonal antibody

Catalog No: #62059

Package Size: #62059-1 50ul #62059-2 100ul



Description

Product Name	TNNT1 rabbit polyclonal antibody
Host Species	Rabbit
Clonality	polyclonal
Purification	Antigen affinity purification
Applications	IHC WB
Species Reactivity	Human Mouse
Immunogen Description	Fusion protein of human TNNT1
Conjugates	Unconjugated
Other Names	troponin T1, slow skeletal type; ANM; TNT; NEM5; STNT; TNTS
Concentration	0.7 mg/ml
Formulation	pH7.4 PBS, 0.05% NaN3, 40% Glycerol
Storage	Store at -20Λ C

Application Details

WB 1:500-2000;IHC 1:25-50

Note: This product is for in vitro research use only and is not intended for use in humans or animals.