

Cubilin Rabbit pAb

Catalog No: #62067

Package Size: #62067-1 50ul #62067-2 100ul


Description

Product Name	Cubilin Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC;IF;WB;ELISA
Species Reactivity	Human;Rat;Mouse
Specificity	Cubilin Polyclonal Antibody detects endogenous levels of Cubilin protein.
Immunogen Description	Synthesized peptide derived from the N-terminal region of human Cubilin.
Target Name	CUBN
Other Names	CUBN ; IFCR ; Cubilin ; 460 kDa receptor ; Intestinal intrinsic factor receptor ; Intrinsic factor-cobalamin receptor ; Intrinsic factor-vitamin B12 receptor
Calculated MW	400kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000; IHC 1:100-1:300; ELISA 1:40000; IF 1:50-200

Images



Immunohistochemical analysis of paraffin-embedded human tonsil.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.