

## DAPL1 Rabbit pAb

Catalog No: #62069

Package Size: #62069-1 50ul #62069-2 100ul


## Description

|                       |                                                                                                                       |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name          | DAPL1 Rabbit pAb                                                                                                      |
| Host Species          | Rabbit                                                                                                                |
| Clonality             | Polyclonal                                                                                                            |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications          | IHC;IF                                                                                                                |
| Species Reactivity    | Human;Mouse                                                                                                           |
| Specificity           | This antibody detects endogenous levels of DAPL1 at Human/Mouse                                                       |
| Immunogen Description | Synthesized peptide derived from human DAPL1                                                                          |
| Target Name           | DAPL1 EEDA                                                                                                            |
| Calculated MW         | 12kD                                                                                                                  |
| Concentration         | 1 mg/ml                                                                                                               |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| Storage               | Store at -20°C/1 year                                                                                                 |

## Application Details

IHC 1:50-200; IF 1:50-200

## Images



Immunohistochemical analysis of paraffin-embedded human uterus.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.