

Band 3/AE1 Rabbit pAb

Catalog No: #62084

Package Size: #62084-1 50ul #62084-2 100ul

Description

Product Name	Band 3/AE1 Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of Band 3/AE1 at Human, Mouse,Rat
Immunogen Description	Synthesized peptide derived from human Band 3/AE1
Target Name	SLC4A1 AE1 DI EPB3
Calculated MW	100kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.