

ULK1 (Phospho Ser555) Rabbit pAb

Catalog No: #62086

Package Size: #62086-1 50ul #62086-2 100ul

Description

Product Name	ULK1 (Phospho Ser555) Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse
Specificity	This antibody detects endogenous levels of Human Mouse ULK1 (phospho-Ser555)
Immunogen Description	Synthesized phospho peptide around human ULK1 (Ser555)
Target Name	ULK1 KIAA0722
Other Names	Serine/threonine-protein kinase ULK1 ; Autophagy-related protein 1 homolog ; ATG1 ; hATG1 ; Unc-51-like kinase 1
Calculated MW	115kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.