

PPAR alpha (phospho-Ser12) Antibody

Catalog No: #62087

Package Size: #62087 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PPAR alpha (phospho-Ser12) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody is from purified rabbit serum by affinity purification via sequential chromatography on phospho-peptide and non-phospho-peptide affinity columns.
Applications	WB;IHC
Species Reactivity	Human; Mouse; Rat
Specificity	The antibody detects endogenous level of PPAR alpha only when phosphorylated at Ser12
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human PPARG, corresponding to amino acid residues around phosphorylated Ser12.
Target Name	MGC2237; MGC2452; NR1C1
Calculated MW	Predicted band size: 52 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:1000-3000;IHC 1:50-1:200

Images



Formalin-fixed, paraffin-embedded Mouse testis tissue for PPAR alpha (phospho-Ser12) Antibody was diluted at 1:100 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.