

Beclin 1 (Phospho-Ser234) Antibody

Catalog No: #62093

Package Size: #62093 100ul

Description

Product Name	Beclin 1 (Phospho-Ser234) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody is from purified rabbit serum by affinity purification via sequential chromatography on phospho-peptide and non-phospho-peptide affinity columns.
Applications	ELISA(peptide)
Species Reactivity	Human, Mouse, Rat
Specificity	The antibody detects endogenous level of Beclin 1 only when phosphorylated at Ser234
Immunogen Type	Peptide
Immunogen Description	Peptide sequence around phosphorylation site of Ser234 (E-Y-S(p)-E-F) derived from Human Beclin 1.
Conjugates	Unconjugated.
Target Name	Beclin 1
Calculated MW	52kD(Calculated)
Concentration	1mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ELISA(peptide) 1:20000-1:40000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.