

## RELM beta antibody

Catalog No: #62101

Package Size: #62101 100ul

## Description

|                       |                                                                                                                       |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name          | RELM beta antibody                                                                                                    |
| Host Species          | Rabbit                                                                                                                |
| Clonality             | Polyclonal                                                                                                            |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications          | IHC                                                                                                                   |
| Species Reactivity    | Human                                                                                                                 |
| Immunogen Type        | Peptide                                                                                                               |
| Immunogen Description | peptide derived from Human RELM beta..                                                                                |
| Target Name           | RELM beta                                                                                                             |
| Calculated MW         | 12kDa                                                                                                                 |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| Storage               | Store at +4°C for short term. Store at -20°C for long term. Avoid freeze/thaw cycle.                                  |

## Application Details

IHC 1:50-1:200

## Background

Probable hormone

Note: This product is for in vitro research use only and is not intended for use in humans or animals.