

## PRRC2C Antibody

Catalog No: #62104

Package Size: #62104 100ul

## Description

|                       |                                                                                 |
|-----------------------|---------------------------------------------------------------------------------|
| Product Name          | PRRC2C Antibody                                                                 |
| Host Species          | Rabbit                                                                          |
| Clonality             | Polyclonal                                                                      |
| Purification          | affinity purified by Protein A                                                  |
| Applications          | IHC;ICC/IF                                                                      |
| Species Reactivity    | (predicted: Human, Mouse, Rat)                                                  |
| Specificity           | The antibody detects endogenous level of total PRRC2C protein                   |
| Immunogen Type        | Peptide                                                                         |
| Immunogen Description | KLH conjugated synthetic peptide derived from human PRRC2C                      |
| Formulation           | 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.               |
| Storage               | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

## Application Details

IHC 1:100-500;ICC/IF 1:100-500

Note: This product is for in vitro research use only and is not intended for use in humans or animals.