

ZNF687 Antibody

Catalog No: #62105

Package Size: #62105 100ul


Description

Product Name	ZNF687 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB;IHC;IF
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total ZNF687 protein
Immunogen Type	Peptide
Immunogen Description	KLH conjugated synthetic peptide derived from human ZNF687
Calculated MW	Predicted band size: 129 kDa
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-2000;IHC 1:100-500;IF 1:100-500

Images



All lanes : ZNF687 Antibody at 1/300 dilution

Lane 1 : A549 whole cell lysates

Lane 2 : U251 whole cell lysates

Lysates/proteins at 20 ug per lane.


Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/5000 dilution

Predicted band size: 129 kDa

Observed band size: 129 kDa

Exposure time: 25 seconds



Formalin-fixed, paraffin-embedded Human kidney stained for ZNF687 Antibody was diluted at 1:400 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.