

Transmembrane protein 30A Antibody

Catalog No: #62109

Package Size: #62109 100ul


Description

Product Name	Transmembrane protein 30A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC;IF
Species Reactivity	Mouse;Rat
Specificity	The antibody detects endogenous level of total Transmembrane protein 30A protein
Immunogen Type	Peptide
Immunogen Description	KLH conjugated synthetic peptide derived from human Transmembrane protein 30A
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.


Application Details

IHC 1:100-500;IF 1:100-500

Images



Formalin-fixed, paraffin-embedded mouse kidney for Transmembrane protein 30A Antibody was diluted at 1:200 dilution in immunohistochemical analysis.



Formalin-fixed, paraffin-embedded rat brain for Transmembrane protein 30A Antibody was diluted at 1:200 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.