

PADI1 Antibody

Catalog No: #62113

Package Size: #62113 100ul

Description

Product Name	PADI1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	(predicted: Human, Mouse, Rat)
Specificity	The antibody detects endogenous level of total PADI1 protein
Immunogen Type	protein
Immunogen Description	Recombinant human PADI1 protein
Calculated MW	Predicted band size: 75 kDa
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.