

EBNA 3B Antibody

Catalog No: #62114

Package Size: #62114 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EBNA 3B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB;IHC;IF
Species Reactivity	(predicted: Epstein Barr Virus)
Specificity	The antibody detects endogenous level of total EBNA 3B protein
Immunogen Type	Peptide
Immunogen Description	KLH conjugated synthetic peptide derived from Epstein Barr Virus EBNA
Calculated MW	Predicted band size: 103 kDa
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-2000;IHC 1:100-500;IF 1:50-200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.