

# Hepatitis B Virus Surface Antigen Antibody

Catalog No: #62116

Package Size: #62116 100ul



## Description

Product Name	Hepatitis B Virus Surface Antigen Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	(predicted: Hepatitis B virus(HBVD2), Hepatitis B virus(HBVC5))
Specificity	The antibody detects endogenous level of total Hepatitis B Virus Surface Antigen protein
Immunogen Type	Peptide
Immunogen Description	KLH conjugated synthetic peptide derived from HBSAG_HBVD2 pre S1 protein
Calculated MW	Predicted band size: 31/44 kDa
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB 1:500-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.