

Staphylococcus aureus Antibody

Catalog No: #62117

Package Size: #62117 100ul

Description

Product Name	Staphylococcus aureus Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ELISA
Species Reactivity	(predicted: Staphylococcus aureus)
Specificity	The antibody detects endogenous level of total Staphylococcus aureus protein
Immunogen Type	protein
Immunogen Description	whole cell protein of Staphylococcus aureus Rosenbach tropina
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ELISA 1:5000-10000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.