

HSV1 + HSV2 gD Antibody

Catalog No: #62119

Package Size: #62119 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HSV1 + HSV2 gD Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC;IF
Species Reactivity	(predicted: HSV1 + HSV2)
Specificity	The antibody detects endogenous level of total HSV1 + HSV2 gD protein
Immunogen Type	Peptide
Immunogen Description	KLH conjugated synthetic peptide derived from HSV1 + HSV2 gD (type 1/strain HZT)
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC 1:100-500;IF 1:100-500

Note: This product is for in vitro research use only and is not intended for use in humans or animals.