

Shigella tropina Antibody

Catalog No: #62123

Package Size: #62123 100ul

Description

Product Name	Shigella tropina Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ELISA
Species Reactivity	Shigella
Immunogen Type	Protein
Immunogen Description	Whole cell protein of Shigella tropina
Conjugates	Unconjugated
Target Name	Shigella tropina
Concentration	1mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ELISA 1:5000-10000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.