

PALM3 Antibody

Catalog No: #62128

Package Size: #62128 100ul

Description

Product Name	PALM3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total PALM3 protein.
Immunogen Type	Peptide
Immunogen Description	KLH conjugated synthetic peptide derived from human PALM3: 1-100/673
Conjugates	Unconjugated.
Target Name	PALM3
Calculated MW	72kD(Calculated)
Concentration	1mg/ml
Formulation	Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.