

NRAMP1 Antibody

Catalog No: #62130

Package Size: #62130 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NRAMP1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC,IF
Species Reactivity	Rat (predicted: Human,Rabbit,Cow)
Specificity	The antibody detects endogenous level of total NRAMP1 protein.
Immunogen Type	Peptide
Immunogen Description	KLH conjugated synthetic peptide derived from human NRAMP1: 371-470/550
Conjugates	Unconjugated.
Target Name	NRAMP1
Calculated MW	60kD(Calculated)
Concentration	1mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC 1:100-500, IF 1:100-500

Images



Formalin-fixed,paraffin-embedded Rat spleen tissue stained for NRAMP1 at 1/100 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.