

## Adenovirus 5 E1A Antibody

Catalog No: #62133

Package Size: #62133 100ul


## Description

|                       |                                                                                   |
|-----------------------|-----------------------------------------------------------------------------------|
| Product Name          | Adenovirus 5 E1A Antibody                                                         |
| Host Species          | Rabbit                                                                            |
| Clonality             | Polyclonal                                                                        |
| Purification          | affinity purified by Protein A                                                    |
| Applications          | WB                                                                                |
| Specificity           | The antibody detects endogenous level of total Adenovirus 5 E1A protein.          |
| Immunogen Type        | Protein                                                                           |
| Immunogen Description | KLH conjugated synthetic peptide derived from human Adenovirus 5 E1A: 181-280/289 |
| Conjugates            | Unconjugated.                                                                     |
| Target Name           | Adenovirus 5 E1A                                                                  |
| Calculated MW         | 32kD(Calculated)                                                                  |
| Concentration         | 1mg/ml                                                                            |
| Formulation           | 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.                 |
| Storage               | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.   |


## Application Details

WB 1:1000-2000

## Images



All lanes : Adenovirus 5 E1A Antibody at 1/1k dilution  
Lane 1 : Recombinant ADV 5E1A protein, His  
Lysates/proteins at 20 ug per lane.  
Secondary  
All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution  
Predicted band size: 32 kDa  
Observed band size: 40 kDa  
Exposure time: 13 seconds



All lanes : Adenovirus 5 E1A Antibody at 1/1k dilution  
Lane 1 : Recombinant Adenovirus 5 E1A protein (50ng)  
Lysates/proteins at 20 ug per lane.  
Secondary  
All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution  
Predicted band size: 32 kDa  
Observed band size: 40 kDa  
Exposure time: 15 seconds

Note: This product is for in vitro research use only and is not intended for use in humans or animals.