

eNOS (Phospho-Ser632) Antibody

Catalog No: #62134

Package Size: #62134 100ul

Description

Product Name	eNOS (Phospho-Ser632) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB,IHC
Species Reactivity	predicted: Human,Mouse,Rat,Rabbit,Pig,Sheep,Dog
Specificity	The antibody detects endogenous level of eNOS only when phosphorylated at Ser632
Immunogen Type	Peptide
Immunogen Description	Peptide sequence around phosphorylation site of Ser632 (K-E-S(p)-S-N) derived from eNOS .
Conjugates	Unconjugated.
Target Name	eNOS
Calculated MW	133kD(Calculated)
Concentration	1mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:1000-2000,IHC 1:100-1:500

Note: This product is for in vitro research use only and is not intended for use in humans or animals.