

Hepatitis B Virus Surface Antigen antibody

Catalog No: #62135

Package Size: #62135 100ul

Description

Product Name	Hepatitis B Virus Surface Antigen antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity purified by Protein A.
Applications	WB, ELISA
Species Reactivity	Hepatitis B virus(HBVD2),Hepatitis B virus(HBVC5)
Immunogen Type	Peptide
Immunogen Description	peptide derived from HBSAG_HBVD2 pre S1 protein: 31-100/389
Target Name	Hepatitis B Virus Surface Antigen
Calculated MW	31/44 kDa
Concentration	1mg/ml
Formulation	Liquid in 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Store at +4°C for short term. Store at -20°C for long term. Avoid freeze/thaw cycle.

Application Details

WB 1:500-2000, ELISA 1:5000-10000

Background

Hepatitis B Virus (HBV) infection induces a disease state characterised by liver damage, inflammation and viral persistence. Infection also increases the risk of hepatocellular carcinoma. HBV belongs to the Hepadnaviridae family of viruses. Its genome consists of partially double stranded circular DNA. The DNA is enclosed in a nucleocapsid, or core antigen (HBcAg), which is surrounded by a spherical envelope (surface antigen or HBsAg). The core antigen shares its sequences with the e antigen (HBeAg) but no cross reactivity between the two proteins has been observed. The HBV genome also encodes a DNA polymerase that also acts as a reverse transcriptase. Hepatitis B infection is normally diagnosed from serological tests that detect HBsAg but as the disease progresses this antigen may no longer be present in the blood and tests for HBcAg are used. If HBsAg can be detected in the blood for longer than six months, chronic hepatitis B is diagnosed. The antigenic determinant of the protein moiety of the HBsAg determines specific characteristics of different serotypes and provides the basis of immunodetection. HBsAg has antigenic heterogeneity, specifically, two pairs of sub specific determinants, d/y and w/r allow the following combinations: adw, ayw, adr, ayr.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.