

Cytomegalovirus pp65 antibody

Catalog No: #62152

Package Size: #62152 100ul


Description

Product Name	Cytomegalovirus pp65 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC, IF
Species Reactivity	Human
Immunogen Description	KLH conjugated synthetic peptide derived from HCMV pp65: 251-350/561
Target Name	HCMV pp65
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P 1:100-500, IHC-F 1:100-500, IF 1:100-500

Images



Immunohistochemical analysis of paraffin-embedded human lung carcinoma with Cytomegalovirus pp65 antibody at 1:200.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.