

HIV2 gp36 antibody

Catalog No: #62154

Package Size: #62154 100ul

Description

Product Name	HIV2 gp36 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB, ELISA
Species Reactivity	HIV2
Immunogen Description	KLH conjugated synthetic peptide derived from HIV2 gp36: 78-125/858
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-2000, ELISA 1:5000-10000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.