

Hepatitis B Virus Surface Antigen Antibody

Catalog No: #62157

Package Size: #62157 100ul

Description

Product Name	Hepatitis B Virus Surface Antigen Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	Immunogen affinity purified
Applications	ELISA
Species Reactivity	Hepatitis B virus
Specificity	The antibody detects endogenous level of total Hepatitis B Virus Surface Antigen protein
Immunogen Description	Recombinant Hepatitis B Virus Surface Antigen
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.