

TRDC Antibody

Catalog No: #62161

Package Size: #62161 100ul

Description

Product Name	TRDC Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB;IF
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total TRDC protein.
Immunogen Type	Peptide
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the Central region of human TRDC. The exact sequence is proprietary.
Conjugates	Unconjugated
Target Name	TRDC
Calculated MW	27kDa
Concentration	1mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200; IF:1:50-1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.