

TRDC Antibody

Catalog No: #62162

Package Size: #62162 100ul


Description

Product Name	TRDC Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total TRDC protein.
Immunogen Type	Peptide
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the Central region of human TRDC. The exact sequence is proprietary.
Conjugates	Unconjugated
Target Name	TRDC
Calculated MW	27kDa
Concentration	1mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-1:1000

Images



All lanes : TRDC Antibody at 1/500 dilution
 Lane 1 : MDAMB453 whole cell lysates
 Lysates/proteins at 20 ug per lane.
 Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
 Predicted band size: 17 kDa Observed band size: 27 kDa
 Exposure time: 3 seconds

Note: This product is for in vitro research use only and is not intended for use in humans or animals.