

TCR alpha/beta Antibody

Catalog No: #62164

Package Size: #62164 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TCR alpha/beta Antibody
Host Species	Mouse
Purification	The antibody was purified by affinity chromatography.
Applications	IF, FC
Species Reactivity	Human
Conjugates	Unconjugated
Target Name	TRAC; TRBC1
Formulation	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, and 0.02% sodium azide.
Storage	Store at 4°C in dark for 6 months

Note: This product is for in vitro research use only and is not intended for use in humans or animals.