

Semaphorin 4B Antibody

Catalog No: #62166

Package Size: #62166 100ul


Description

Product Name	Semaphorin 4B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total Semaphorin 4B protein.
Immunogen Type	Peptide
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the C-terminal region of human Semaphorin 4B. The exact sequence is proprietary.
Conjugates	Unconjugated
Target Name	SEMA4B
Calculated MW	105kDa
Concentration	1mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-1:1000

Images



All lanes : Semaphorin 4B Antibody at 1/1K dilution
 Lane 1 : A549 whole cell lysates
 Lane 2 : DU145 whole cell lysates
 Lysates/proteins at 20 ug per lane.
 Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
 Predicted band size: 92 kDa Observed band size: 105 kDa
 Exposure time: 12 seconds

Note: This product is for in vitro research use only and is not intended for use in humans or animals.