

DLG5 Antibody

Catalog No: #62180

Package Size: #62180 100ul



Description

| | |
|-----------------------|---|
| Product Name | DLG5 Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications | WB |
| Species Reactivity | Human |
| Specificity | The antibody detects endogenous level of total DLG5 protein |
| Immunogen Type | Peptide |
| Immunogen Description | Synthesized peptide derived from human DLG5 |
| Target Name | DLG5; KIAA0583; PDLG; |
| Calculated MW | Predicted band size: 211 kDa |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

WB 1:500-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.