

TMEM24 Antibody

Catalog No: #62196

Package Size: #62196 100ul

Description

Product Name	TMEM24 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human; Mouse
Specificity	The antibody detects endogenous level of total TMEM24 protein
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from human TMEM24
Target Name	C2CD2L; KIAA0285; TMEM24; DLNB23;
Calculated MW	Predicted band size: 78 kDa
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.