

HECTD3 Antibody

Catalog No: #62200

Package Size: #62200 100ul

Description

Product Name	HECTD3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human; Mouse
Specificity	The antibody detects endogenous level of total HECTD3 protein
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from human HECTD3
Target Name	HECTD3;
Calculated MW	Predicted band size: 95 kDa
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.