

## APM2 Antibody

Catalog No: #62206

Package Size: #62206 100ul


## Description

|                       |                                                                                                                       |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name          | APM2 Antibody                                                                                                         |
| Host Species          | Rabbit                                                                                                                |
| Clonality             | Polyclonal                                                                                                            |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications          | IHC;IF                                                                                                                |
| Species Reactivity    | Human                                                                                                                 |
| Specificity           | The antibody detects endogenous level of total APM2 protein                                                           |
| Immunogen Type        | Peptide                                                                                                               |
| Immunogen Description | Synthesized peptide derived from human APM2                                                                           |
| Target Name           | APM2; C10orf116                                                                                                       |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| Storage               | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.                                       |

## Application Details

IHC 1:50-200; IF 1:50-200

## Images



Formalin-fixed, paraffin-embedded human tonsil for APM2  
Antibody was diluted at 1:200 dilution in  
immunohistochemical analysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.